

flexmobil comfort



 **EUROFLEX**



flexmobil^{comfort} - stable Comfort

Flexmobil Comfort is a power wheelchair for indoor use. Very easy to drive thanks to the driving wheels placement under the seat in the center of the chassis, which makes it possible to rotate the wheelchair around its own center. Together with the 4 swivel wheels this makes the wheelchair very stable at all times. Flexmobil Comfort is equipped with the Comfort seat unit, developed for users who spend a lot of time in their wheelchair.

The dependable and popular Flexmobil can be ordered with this comfortable seat unit Comfort.

A lot of customers appreciate the simple and reliable construction of Flexmobil, together with the easy handling and stable construction. The already good comfort has reached an even higher level with the Comfort seat unit.



In the Comfort seat you can sit a long time without fatigue.



The chassis with 6 wheels chassis and driving on the center wheels makes Flexmobil very easy to drive with minimal driving space..



The quick release of the back support makes the wheel chair easy to transport in the trunk of a vehicle, e.g. station wagon.



power seatlift



power tilt



power back



power legsupport

Standard features:

- chassis - flexmobil - center wheel drive
- wheels - drive wheel and castors - filled
- comfort seat - 45 cm, seat depth 40-61 cm
- seat lift - power
- back recline - manual - gas assisted
- seat tilt - manual - screw (onetime adjust)
- arm support - soft - fold back
- elektronics - Dynamic - G90T

Accessories: (a selection)

- comfort seat unit - 40 cm, 53 cm -low back-51cm
- comfort seat unit - 45 cm - high back - 64 cm
- seat tilt - power
- back angle - power
- leg support - single foot plate
- leg support - divided manual or power
- head rest
- trunk support
- side support
- arm support - extra soft
- abductor
- table, tray

Technical data:

Class	A
Obstacles	max 30 mm
Driving	center wheels
Seat, width Comfortsits	40, 45, 53 cm
Back, height Comfortsits	51 low, 64 high cm
Seat, depth -from back to seat front	40-61 cm, adjustable
Seat, width- between arm support	46-56 cm, adjustable
Seat, height - to seat top	45-63 cm, alt. 50-75 cm
Sits, tilt, power	-4° forward to +16° backward
Back, recline	4° forward to -30° backward
Fabric	leather
Armrest, height	10-30 cm
Armrest, size	31x8 cm
Total lenght - incl. footrest	106 cm
Total lenght - chassis	82 cm
Total width - chassis	58 cm
Total weight	85 kg
Max user weight	125 kg
Turning space	106 cm
Max driving time	8,5 h
Max speed	4 km/h
Electronics	Dynamic G90T

