



**VOLARIS**  
SWEDEN

# Volaris S7

## Light and flexible

The Volaris S7 is a newly developed light and flexible rollator with functions for everybody. The collapsible X-type construction facilitates the fold up and makes the rollator stand steady in folded position. This means that the rollator requires less space and is easy to transport. As both the seat and the handles are easily vertically adjustable, the sole version of the rollator is suitable for varying groups of users. The rollator is equipped with VABS brake for efficient braking. The completely new and unique parking brake in form of a loop is plainly and user-friendly. Volaris S7 is manufactured in Sweden from re-usable material. Design and construction patent applied for. The rollator has been developed in accordance with the applicable standard for walking aids, ISO 11199-2:2006.

### Technical specification:

Weight: 7,2 kg

Weight capacity: 140kg



Mobility for everyone™

[www.volaris.se](http://www.volaris.se)

Artnr# 88 104 02 | 091113

