

Low Seat height

Volaris S7 Low is one of the shortest rollators on the market. The innovative features of the rollator include a new and unique parking brake in the form of a user friendly loop. The Volaris Advanced Braking System is also distinctive in that the brakes can be applied with minimal effort. Other features include the flexible 'comfort frame' which allow the rollator to remain stable on uneven surfaces.

The frame flexes enabling the wheels to remain on the ground when going over bumps. The Xfolded design facilitates folding and allows the rollator to stand steady in the folded position. The light weight enables you to lift the walker into the back of the car.

Accessories can be fitted to suit individual requirements.

Volaris S7 Low is made in Sweden of recyclable material. The rollator is developed according to the ISO standard ISO 11199-2:2005.

Technical specifications:

Weight: 7,1 kg

Personal weight capacity: 140kg



Mobility for everyone™

www.volaris.se

090810

