


IMPORTANT SAFETY NOTICE PIXI SULKY STROLLER

1 June 2015

Dear Pixi Sulky Stroller user.

We have noticed that there is a small risk that the handle joint can break, if the instructions we provide in the manual are not observed. This is a clarification of how to handle the Stroller safely when lifting or on stairs. The new manual has been updated in accordance with this instruction.

- **NOTE!** When using on stairs and when lifting, the handle must be folded down. Grip the chassis inside/on the handle joint for safe handling of stroller and child.



This notice must be forwarded to all those concerned within your organisation, or to another organisation that has received this product that could constitute a risk. If you have received this notice at home, it has been sent to you in close cooperation with your health care provider. Your personal details are processed confidentially and will not be used for any other purpose, nor will they be stored by Eurovema.

The undersigned hereby declares that this notice has been reported to the appropriate regulatory authority.

Stefan Helgebring
Technical Manager
Eurovema AB
Baldersvägen 38
SE-332 35 Gislaved
Phone +46 371 390 100



Eurovema AB
Baldersvägen 38
SE-332 35 Gislaved
Sweden

Telephone:
+46 (0)371-39 01 00
Telefax:
+46 (0)371-189 82

E-mail:
info@eurovema.se
Internet:
www.eurovema.se

Bankgiro:
5008-7261
Danske Bank

Org.No:
556592-4411
VAT.No:
SE556592441101



VOLARIS
SWEDEN