

Order

Quantity

Quotation

Invoice adress

Delivery adress

Tel:  Fax:

Order nr:  Date:

Square =standard

E-mail:

## Model

### Chairs with gas assisted seat lift

- |   |          |
|---|----------|
| <input type="checkbox"/> Manual seat angle, low back rest           | 241-XX*  |
| <input type="checkbox"/> Manual seat angle, high back rest          | 242-XX*  |
| <input type="checkbox"/> Coxit sea, man. seat angle, low back rest  | 241C-XX* |
| <input type="checkbox"/> Coxit seat, man. seat angle, high backrest | 242C-XX* |

Low backrest=35x30cm  
High backrest=35x44 cm

\*XX is replaced with **Art.nr.** for color, se below.  
JEA is electric seat angle

### Chairs with electric seat lift

- |  |              |
|--|--------------|
| <input type="checkbox"/> Manual seat angle, low back rest            | 291-XX*      |
| <input type="checkbox"/> Manual seat angle, high back rest           | 292-XX*      |
| <input type="checkbox"/> Electric seat angle, low back rest          | 291-JEA-XX*  |
| <input type="checkbox"/> Electric seat angle, high back rest         | 292-JEA-XX*  |
| <input type="checkbox"/> Electric tilt, low back rest                | 291-2JEA-XX* |
| <input type="checkbox"/> Electric tilt, high back rest               | 292-2JEA-XX* |
| <input type="checkbox"/> Coxit seat, man. seat angle low back rest   | 291C-XX*     |
| <input type="checkbox"/> Coxit seat, man. Seat angle, high back rest | 292C-XX*     |

### Color (upholstery)

- |                                     |    |
|-------------------------------------|----|
| <input type="checkbox"/> Black dot  | 33 |
| <input type="checkbox"/> Gray video | 11 |
| <input type="checkbox"/> Blue video | 12 |

### Chassis

- |  |      |
|--|------|
| <input type="checkbox"/> Standard 48x53 cm                 | -    |
| <input type="checkbox"/> Flexbase 54x55 cm                 | 2401 |
| <input type="checkbox"/> H-chassis 50x60cm (2-wheel brake) | 2402 |
| <input type="checkbox"/> H-chassis 50x60cm (4-wheel brake) | 2403 |

### Wheels

- |                                 |      |
|---------------------------------|------|
| <input type="checkbox"/> 100 mm | -    |
| <input type="checkbox"/> 125 mm | 2404 |
| <input type="checkbox"/> 150mm  | 2504 |

### Electric brakes

- |   |        |
|---|--------|
| <input type="checkbox"/> Standard chassis | 2909   |
| <input type="checkbox"/> Flexbase chassis | 2909-1 |
| <input type="checkbox"/> H-chassis        | 2910   |

### Seat hight and user weight (gas assisted seat lift) Art.nr.

- |  |       |
|--|-------|
| <input type="checkbox"/> 43-57 cm, 45-140 kg                           | -     |
| <input type="checkbox"/> 43-57 cm, 90-140 kg                           | 2431  |
| <input type="checkbox"/> 45-65 cm, 45-140 kg, long                     | 2433  |
| <input type="checkbox"/> 38-46 cm, 45-140 kg, short                    | 2432  |
| <input type="checkbox"/> Seat rotation                                 | 2406  |
| <input type="checkbox"/> Seat hight adjustment under right arm support | 2433H |
| <input type="checkbox"/> Seat hight adjustment under left arm support  | 2433V |

The position of the gas spring in the chassis is adjustable

### Seat hight and user weight (electric seat lift) Art.nr.

- |   |      |
|---|------|
| <input type="checkbox"/> 45-70 cm, 140 kg | -    |
| <input type="checkbox"/> 40-58 cm, 140 kg | 2906 |
| <input type="checkbox"/> Seat rotation    | 2907 |

125kg with electric seat angle or tilt

The position of the actuator in the chassis is adjustable

<b>Seat plate</b>	<b>Art.nr.</b>	
<input type="checkbox"/> 46x46 cm	-	
<input type="checkbox"/> 43x43 cm	2407-XX	
<input type="checkbox"/> 40x40cm	2408-XX	
<input type="checkbox"/> 40x46 cm	2409-XX	
<input type="checkbox"/> 49x51cm	2410-XX	Seat sizes 43x43 and 46x46 cm comes with removable cover. Coxit seat available only in size 46x46 and 49x51 cm.
<input type="checkbox"/> Plane seat	2412	

<b>Armrest LEFT</b>	<b>Art.nr.</b>	<b>Armrest RIGHT</b>	<b>Art.nr.</b>
<input type="checkbox"/> Polyurethane 30 cm with plush	-	<input type="checkbox"/> Polyurethane 30 cm with plush	-
<input type="checkbox"/> Soft 30 cm	2417V	<input type="checkbox"/> Soft 30 cm	2417H
<input type="checkbox"/> Soft 40 cm	2418V	<input type="checkbox"/> Soft 40 cm	2418H
<input type="checkbox"/> Hemiplegic armrest	2422V	<input type="checkbox"/> Hemiplegic armrest	2422H
<input type="checkbox"/> Comfort 30 cm	2430V	<input type="checkbox"/> Comfort 30 cm	2430H
<input type="checkbox"/> Foldeble	2419V	<input type="checkbox"/> Foldeble	2419H

<b>Leg support LEFT</b>	<b>Art.nr.</b>	<b>Leg support RIGHT</b>	<b>Art.nr.</b>
<input type="checkbox"/> Devided leg support,	224001-2000V	<input type="checkbox"/> Devided leg support,	224001-2000H
<input type="checkbox"/> Amputation support	224001-4100V	<input type="checkbox"/> Amputation support	224001-4100H
<input type="checkbox"/> Foot plate	224001-2019		
<input type="checkbox"/> Foot plate universal	259001-2004		
<input type="checkbox"/> Footring foldeble	624001-2020		
<input type="checkbox"/> Footring foldeble, Flexbase	684001-2020-xl		
<b>Adapter for other manufacturers leg supports</b>		<b>Adapter for other manufacturers leg supports</b>	
<input type="checkbox"/> ETAC	624001-2104V	<input type="checkbox"/> ETAC	624001-2104H
<input type="checkbox"/> Combi hold	224001-2002	<b>Must be included when leg support is choosed</b>	

<b>Seat cover (inc. protection)</b>	<b>Art.nr.</b>	<b>Seat cover (soft)</b>	
<input type="checkbox"/> 46x46 cm	6243xx-3051	<input type="checkbox"/> 46x46 cm	6242xx-3051-SR3
<input type="checkbox"/> 43x43 cm	6240xx-3051	<input type="checkbox"/> 43x43 cm	6240xx-3051-SR3
<input type="checkbox"/> 40x40 cm	6241xx-3051	<input type="checkbox"/> 40x40 cm	6241xx-3051-SR3
<input type="checkbox"/> 40x46 cm	6242xx-3051	<input type="checkbox"/> 40x46 cm	6242xx-3151-SR3
<input type="checkbox"/> 49x51 cm	6244xx-3051	<input type="checkbox"/> 49x51 cm	6244xx-3051-SR3
<input type="checkbox"/> Coxit	6243xx-3051C	<input type="checkbox"/> Coxit	6243xx-3051C-SR
<b>Back rest cover (inc. protection)</b>		<b>Back rest cover (inc. protection)</b>	
<input type="checkbox"/> 35x30 cm	6241xx-6151	<input type="checkbox"/> 35x44cm	6242xx-6151

**Accessories****Art.nr.****Art.nr.**

<input type="checkbox"/> Neck support Comfort (only with high back rest)	299001-4700	<input type="checkbox"/> Side support 10x24 cm LEFT	224047-4500V
<input type="checkbox"/> Neck support U-shape (only with high back rest)	234001-4700	<input type="checkbox"/> Side support 10x24 cm RIGHT	224047-4500H
<input type="checkbox"/> Table, foldble sideways	224001-4900	<input type="checkbox"/> Side support 10x32 cm LEFT	224047-4501V
<input type="checkbox"/> Pelvic seat belt	624047-4400	<input type="checkbox"/> Side support 10x32 cm RIGHT	224047-4501H
<input type="checkbox"/> 4-point seat belt	624001-4100	<input type="checkbox"/> Abductor, in black Skai	224047-4800
<input type="checkbox"/> Trunk support LEFT	224047-4600V	<input type="checkbox"/> Push handle, adjustable in hight and angle	224001-2050
<input type="checkbox"/> Trunk support RIGHT	224047-4600H	<input type="checkbox"/> Brake lever extra long 50cm	2423
<input type="checkbox"/> Back rest, low, extra wide, 40x30 cm	2414-xx	<input type="checkbox"/> Brake lever center	2427
<input type="checkbox"/> High back rest, extra wide, 40x44 cm	2415-xx	<input type="checkbox"/> Brake lever foot	2426
<input type="checkbox"/> Back rest, adjustable, Slingback	2416-xx		

**Location of handles****Seat lift handle**

- Right side  
 Left side

**Brake handle**

- Right side  
 Left side

Manual seat angle handle always on left side

**Others**
