

Order

Quantity

Quotation

Invoice adress

Delivery adress

Tel:  Fax:

Order nr:  Date:

E-mail:

**Model**

- 291EXL, seat 60x45 cm, back rest 55x35 cm
- 292EXL, seat 70x50 cm, back rest 65x35 cm

**Art.nr.**

- 391-XX\*
- 392-XX\*

**Standard:**

- 260 kg user weight
- electric seat lift, twin actuators
- electric 4-wheel brake
- Slowrecovery seat
- Adjustable seat angle, +3 to -3 degrees,

**Color** (upholstery)

- Black dot
- Gray video
- Blue video

**Art.nr.**

- 33
- 11
- 12

**Seat hight**

- 48x66 cm, short actuators
- 53x78 cm, long actuators

**Art.nr.**

- 12
- 12

The actuators location in the chassie can be adjusted in 7 steps up to + 7 cm in hight.

**Others**