

## Orderform workchair Forma/SitRite

Square = standard

Order

Quantity

Quotation

Invoice address

Delivery address

A) Undercarriage	No	B) Seat lift	Seat height	No	C) Seat	No
<input type="checkbox"/> Square 100 mm wheels Brake mounted on seat	0	<input type="checkbox"/> Gas piston 350N, <90kg	43-57 cm	0	<input type="checkbox"/> 45x45 cm	0
<input type="checkbox"/> Square 125 mm wheels Brake mounted on seat	1	<input type="checkbox"/> Gas piston 220N, <50kg	43-57 cm	1	<input type="checkbox"/> 40x40 cm	1
<input type="checkbox"/> 4-wheel brake, 100 mm wheels Brake mounted on seat	2	<input type="checkbox"/> Gas piston 450N, <140kg	43-57 cm	2	<input type="checkbox"/> 40x45 cm	2
<input type="checkbox"/> 4-wheel brake, 125 mm wheels Brake mounted on seat	3	<input type="checkbox"/> Gas piston 350N, tall <90kg	45-65 cm	3	<input type="checkbox"/> 50x50 cm	3
<input type="checkbox"/> Square 100 mm wheels Chassis mounted brake	4	<input type="checkbox"/> Gas piston 350N, short <90k	41-49 cm	6	<input type="checkbox"/> Coxit 45x45 cm	4
<input type="checkbox"/> Square 125 mm wheels Chassis mounted brake	5	<input type="checkbox"/> Electric actuator, 250 mm <'	51-76 cm	4	<input type="checkbox"/> SitRite 40x40 cm*	5
<input type="checkbox"/> Electric brake chassis 100 mm	6	<input type="checkbox"/> Electric actuator, 180 mm <'	45-63 cm	5	<input type="checkbox"/> SitRite 40x45 cm*	6
<input type="checkbox"/> Electric brake chassis 125 mm	7	The position of the seat lift can be adjusted in the base which affect high and low seat height.			<input type="checkbox"/> SitRite 45x45 cm*	7
					<input type="checkbox"/> SitRite 45x50 cm*	8
					<input type="checkbox"/> SitRite 50x50 cm*	9
					<input type="checkbox"/> Plane seat	

D) Seat angle	No	E) Back	No	F) Upholstery	No
<input type="checkbox"/> Manual	0	<input type="checkbox"/> Bas-back	0	<input type="checkbox"/> Black*	0
<input type="checkbox"/> Electric-tilt forward-back*	2	<input type="checkbox"/> Bas-back+aktiv back	1	<input type="checkbox"/> Blue*	1
*Tilt: Forward 20° Back 10°		<input type="checkbox"/> Bas-back+komfort back	2	*SitRite in black Atlantic	
<input type="checkbox"/> Electric back	3	<input type="checkbox"/> SitRite 37x43 cm	3		
<input type="checkbox"/> Elektrik-angle forward, electric back	4	<input type="checkbox"/> SitRite 42x45 cm	4		
On next page you will find neck support as accessories.		<input type="checkbox"/> SitRite 47x47 cm	5		
		<input type="checkbox"/> SitRite with lumbar support 37x43cm	6		
		<input type="checkbox"/> SitRite with lumbar support 42x45 cm	7		
		<input type="checkbox"/> SitRite with lumbar support 47x47 cm	8		

G) Armrest	No
<input type="checkbox"/> Standard, foldeble	0
<input type="checkbox"/> No armrest	1
<input type="checkbox"/> Classic pur	2

If other armrest or combinations other than the default is desired, select "no armrest" choose instead armrest on the accessories on the next page.

Fill in the number of the selected design in the boxes, which forms the article number for the chair.

8

A	B	C	D	E	F	G	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

## Location of brake

Right

Left

## Brake system front or rear

Rearbrake

Frontbrake  
Valid only for seat mounted brake.

## Accessories

No accessories

Selected the chair with selected accessories

## Others





## Accessories

### Leg support

	Artno
<input type="checkbox"/> 1/4 footing	800003-1000
<input type="checkbox"/> Footplate, lot	224001-2019
<input type="checkbox"/> Footrest, shared, right	224001-2000H
<input type="checkbox"/> Footrest, shared, left	224001-2000V
<input type="checkbox"/> Footrest bracket for flat seat	224001-2210

## Artno

### Arm rest

	Artno
<input type="checkbox"/> Foldable, soft, black, right	231010-4001H
<input type="checkbox"/> Foldable, soft, black, left	231010-4001V
<input type="checkbox"/> Classic pur, right	224001-4004H
<input type="checkbox"/> Classic pur, left	224001-4004V

## Artno

You need footrest bracket for the footplate, divided footrest and abductor.

## Seat cover (inc. Protection)

	Black	Blue	SitRite
40x40 cm	<input type="checkbox"/> 800401-4040	<input type="checkbox"/> 800402-4040	<input type="checkbox"/> 624155-3051
40x45 cm	<input type="checkbox"/> 800401-4045	<input type="checkbox"/> 800402-4045	<input type="checkbox"/> 624255-3051
45x45 cm	<input type="checkbox"/> 800401-4545	<input type="checkbox"/> 800402-4545	<input type="checkbox"/> 624355-3051
45x50 cm, only for SitRite			<input type="checkbox"/> 624555-3051
50x50 cm	<input type="checkbox"/> 800401-5050	<input type="checkbox"/> 800402-5050	<input type="checkbox"/> 624455-3051
Coxit	<input type="checkbox"/> 800411-4545	<input type="checkbox"/> 800412-4545	

## Headrest

	Black	Blue
Aktiv back	<input type="checkbox"/> 804000-2013	<input type="checkbox"/> 804000-2014
Komfort back	<input type="checkbox"/> 804000-2011	<input type="checkbox"/> 804000-2012
SitRite	<input type="checkbox"/> 804000-2010	

## Back cover

	Black	Blue	SitRite
Bas	<input type="checkbox"/> 801000-4001	<input type="checkbox"/> 801000-4002	37x43 cm <input type="checkbox"/> 624155-6151
Aktiv	<input type="checkbox"/> 802000-4001	<input type="checkbox"/> 802000-4002	42x45 cm <input type="checkbox"/> 624355-6151
Komfort	<input type="checkbox"/> 803000-4001	<input type="checkbox"/> 803000-4002	47x47 cm <input type="checkbox"/> 624255-6151
Headrest covers	<input type="checkbox"/> 804000-4001	<input type="checkbox"/> 804000-4002	

<b>Others</b>	<b>Artno</b>	<b>Others</b>	<b>Artno</b>
<input type="checkbox"/> Hipbelt	224047-4400	<input type="checkbox"/> Seat rotation, gas Not cable braked chassis	224101-3000
<input type="checkbox"/> Abductor	224070-4800	<input type="checkbox"/> Seat rotation, electric Not cable braked chassis	249101-9000
<input type="checkbox"/> Table, foldable, Forma Only for Formas armrest	810000-3020	<input type="checkbox"/> Electric brake, chassismounted brake	800050-2000
<input type="checkbox"/> Table, foldable Classic Only for Classic armrest	224001-4900		
<input type="checkbox"/> Trunk support, black skai, left In combination with comfortback	809000-2011		
<input type="checkbox"/> Trunk support, black skai, right In combination with comfortback	809000-2021		
<input type="checkbox"/> Thigh support, 10x24 cm, black skai, right	800100-1001		
<input type="checkbox"/> Thigh support, 10x24 cm, black skai, left	800100-2001		
<input type="checkbox"/> Thigh support, 10x32 cm, black skai, right	800100-1002		
<input type="checkbox"/> Thigh support, 10x32 cm, black skai, left	800100-2002		
<input type="checkbox"/> Extra trunksupport plate for SitRite (par)	624001-4700		
<input type="checkbox"/> Trunksupport ABC/Sitrite, right	274001-4500H		
<input type="checkbox"/> Trunksupport ABC/Sitrite, left	274001-4500V		