

sprängskiss/exploded view




**Eurovema**  
Creates mobility and ergonomics.

art.nr. SPR-ITO rev: 3/2013-06


**Eurovema**

Creates mobility and ergonomics.






Pos	art.nr	Benämning	Description	Not/Notes
1	60108500	Ito, dyna, latte, small	Ito, cushion, latte, small	komplett/complete
1a	60108600	Ito, dyna, latte, large	Ito, cushion, latte, large	komplett/complete
2	680001-6801	Ito, överdrag, kudde, small	Ito, covers, pilooow, small	
2a	680002-6801	Ito, överdrag, kudde, large	Ito, covers, pilooow, large	
3	70101308	Ito, sidostöd, small	Ito, sidesupport, small	komplett/complete
3	70101309	Ito, sidostöd, large	Ito, sidesupport, large	komplett/complete
4	70101306	Ito, knäled, small	Ito, kneejoint, small	exkl. fotplatta/excl footrest
4a	70101307	Ito, knäled, large	Ito, kneejoint, large	exkl. fotplatta/excl footrest
5	70101304	Ito, fotplatta, small	Ito, footrest, plate, small	komplett/complete
5a	70101305	Ito, fotplatta, large	Ito, footrest, plate, large	komplett/complete
6	70101312	Ito, sits, small	Ito, seat, small	
6a	70101313	Ito, sits, large	Ito, seat, large	



Pos	art.nr	Benämning	Description	Not/Notes
1	70101302	Ito, bakhjul	Ito, rear wheel	komplett/complete
2	70101301	Ito, framhjul	Ito, front wheel	inkl skruv/incl screw
3	70101303	Ito, framhjul med gaffel	Ito, front wheel fork	komplett/complete
4	685001-1400	Ito, bromsfjäder	Ito, brake spring	
5	70101205	Ito, broms, small	Ito, brake, small	komplett/complete
5a	70101206	Ito, broms, large	Ito, brake, large	komplett/complete
6	70101314	Ito, chassi, small	Ito, frame, small	komplett/complete
6a	70101315	Ito, chassi, large	Ito, frame, large	komplett/complete
7	70101210	Ito, tilt, small	Ito, tilt, small	komplett/complete
7a	70101211	Ito, tilt, large	Ito, tilt, large	komplett/complete
8	680001-9999	Ito, handtagsled	Ito, handle part	
9	380001-3100	Ito, handtag, small	Ito, handle, small	
9a	380002-3100	Ito, handtag, large	Ito, handle, large	
10	715001-1210	Ito, bromsgummi	Ito, brake rubber	

# Eurovema

Creates mobility and ergonomics.



Eurovema AB

Baldersvägen 38, SE-332 35 Gislaved, Sweden  
Tel./Phone +46 371 390 100 | Fax. +46 371 189 82  
info@eurovema.se | www.eurovema.se

art.nr. SPR-ITO rev: 3/2013-06