

## Light and flexible



### Kid

Volaris S7 kid are the latest products in the Volaris family. Light and flexible rollators with functions to meet the kids special demand. With its 50 mm narrower and 100 mm shorter frame together with the smaller wheels makes it easy and maneuver in tight spaces and easier to get around.



With frame construction Comfort Frame™ means that the walker takes the soft and come over no-flat surfaces. Rollators folds easily and repository standing with his x-folding design.

The seat and the handle can be easily adjusted in height allows it to accommodate several children with different needs. The walker is equipped with Smart handle. This means that there are no hanging wires that can caught into things, but everything is encapsulated in the frame. Both the process - and the parking brake are ergonomically designed and easily accessible and contribute to high safety.

Volaris S7 Kid is made in Sweden of recyclable materials. Design and construction is patent pending. The walker is designed for existing standards for walking aids, ISO 11199-2:2005

Large accessory program to meet all special needs.



Kid	6,5 kg	53 cm	62 cm	58-70 cm	42 cm	36 cm	150x32 mm	27 cm	60 kg
-----	--------	-------	-------	----------	-------	-------	-----------	-------	-------

## Eurovema

Creates mobility and ergonomics.

Baldersvägen 38, SE-332 35 Gislaved, Sweden  
 Phone. +46 371 390 100 | Fax. +46 371 189 82  
 info@eurovema.se | www.eurovema.se



Rev. 201402|2.