

Making Your walk outdoors a breeze.

Volaris S7 PATROL is a stable rollator with four large wheels made from soft shock absorbing material, making your walk outdoors a breeze. The X-type design facilitates folding and ensures that the rollator stands steadily in the folded position.

The seat and handle height on the Volaris S7 PATROL can be easily adjusted. This means that one walker can be adapted to suit tall or short people. It would be particularly useful in situations where many different users depend on one walker to be adapted to their individual needs, for example in rehabilitation centres.

The backward pointing handles can also be adjusted to suit individual needs with a choice of anatomical as well as cylindrical hand grips. Another unique feature of the S7 PATROL is a patented link system brake function. The travel and parking brakes are very accessible and contribute to a safe and effective braking action.

There are a range of different accessories available to ensure individual user needs are met.

The Volaris S7 PATROL is manufactured in Sweden from reusable material. Patents have been applied for covering both construction and design. The rollator is developed in accordance with the prevailing standard for walking aids, ISO 11199-2:2005

Technical specifications:

Weight: 8,4 kg

Personal weight capacity: 140kg



Mobility for everyone™

www.volaris.se

Artnr# 88 106 02 | 091002

