

*it*



**ALVEMA**





Ergonomics, therapists, designers and technicians have joined forces to produce a pushchair with both the child and carer in mind. The Ito from Alvema which combines simplicity and flexibility, is lightweight and has a contemporary design. The chair has been designed for easy daily use, without making any concession for sitting posture or comfort. Furthermore the Ito can be adapted to meet individual needs, thanks to a well thoughtout selection of accessories.



Ito is easily folded in three steps to a small unit that fits in the car boot.

Accessories (a selection)



shopping basket



swivel or fixed wheels



sun canopy



thermal bag



grip bar



hand brake



abductor



wind and rain cape

Standard features



hip support



seat depth adjust



seat tilt adjust



push handle adjust



Ito front view



Ito side view

Standard features:

- hip support, adjustable
- seat depth, adjustable
- foot rest, adjustable
- push handle, adjustable
- seat covers, soft, removable, washable
- basket
- back extension
- crash hook

Accessories:

- abductor
- grip bar
- sun canope
- thermal bag
- wind and rain cape
- hand brake
- shopping basket
- swivel or fixed front wheels
- head rest
- trunk supports
- vest harness
- tie down brackets

Technical facts:

Size	small	large
Seat, width	20-28 cm	26-34 cm
Seat, depth	22-28 cm	28-35 cm
Back, height	50 / 62 cm	58 / 72 cm
Foot rest, height	14-28 cm	20-36 cm
User weight, maximum	25 kg	40 kg
User height, maximum	120 cm	145 cm
Width, total	58,5 cm	63,5 cm
Size, folded	79,5x58,5x37 cm	84x64,5x37 cm
Weight	12 kg	13 kg
Wheel	250/175 mm	250/175 mm

Material

During manufacturing only non-toxic paint and powder coating is used. All of the upholstery is non-flammable According to standard ISO EN1021-1 and 1021-2. These materials have also been tested for biodegrad ability, standard ISO 10993-5. All seat covers are washable. Ito is crash tested and approved according to standard

Test

Eurovema AB reserves the right to make technical changes.

