

# miniflex<sup>forma</sup>



 **EUROFLEX**



# miniflex<sup>forma</sup> - comfort and ergonomics in mini size

Miniflex Forma is a power wheelchair for indoor use with minimal size and equipped with an advanced workchair seat unit, developed together with experts in ergonomics and seating. The dependable and popular Miniflex can now be ordered with the new and unique seat unit Forma, taken from our new workchair with the same name (Forma). A lot of customers appreciate the simple and reliable construction of Miniflex, together with the easy handling and small size. The already good comfort has now reached an even higher level with the new Forma seat unit.

Miniflex is front wheel driven and the larger front wheels see to that you can pass thresholds easy and the smallness of Miniflex makes you come close to tables and work places. The joint fixed rear wheels keeps all four wheels in contact with the ground all the time.



As standard Miniflex Forma is delivered with the Basic back. Designed to give minimal pressure on the spine and maximum support for the hip.



The Comfort back can be chosen as an accessory. Gives greater support and can be shape to the curve of the users back.



The Active back gives, you good support and still keep the mobility of the users upper part of the torso.



Miniflex Forma Basic in black.



Miniflex Forma Basic in blue.

## Standard:

- miniflex chassi
- forma seat unit
- swing away armsupport
- seat size 45x45 cm - black
- seat angle - manual
- back - model standard - black
- controlbox - Dynamic Shark
- electric actuator - 250 mm stroke

## Accessories: (a selection)

- Comfort Backrest
- Active Backrest
- neck support -in combination with Comfort back
- trunk support
- side suport
- arm support - soft
- abductor
- seat sizes 40x40, 40x45, 50x50, coxit 45x45 cm
- covers - seat and back - inc.protected
- seat tilt - electric
- table
- foot plate
- leg support - separate-manual or electric

## Tecnical facts:

Class	A (indoor use)
Obstacle climb	max 3 cm
Drive	front wheel drive
Seat, width- between armsupports	38-52 cm
Seat, depth -from back to seat front	32-58 cm
Seat, height - to seat top	55-80 cm, 50-68 cm short actuator *
Seat, sizes, WxD , cm	45x45,40x40,40x45, 50x50, 45x45 cm
Seat, angle, manual, forward	0° to +20°
Seat, tilt, EL- forward-backward	-20° to +10°
Back, size, std WxH	37x23 cm
Back, height from seat, std back	21-36 cm
Neck support, size	25x8 cm
Fabric, colors	black, blue
Armsupport, height	19-25 cm
Armsupport, size	31x8 cm
Total length - chassis	74 cm
Total -width - chassis	58 cm
Total weight	60 kg
Max user weight	140 kg
Turning radius	120 cm
Max driving time	8,5 tim
Max speed	4 km/h
Electronics	Dynamic

\* the actuators position in the chassis is adjustable in a screw joint with max. +(-) 5 cm.